


WALKING ON THE BRIGHTER SIDE OF ULTRASOUND IMAGING

MyLab™X
Beyond flexibility








# MyLab™X6 Ultrasound imaging beyond flexibility

**Esaote**'s new **MyLab™X6** makes your workflow so efficient and smooth, it does not only increase productivity, but it also **empowers your clinical performances**.

Take advantage of the 21.5" HD IPS technology LED monitor to get **outstanding image quality** and an unparalleled degree of details from your scans.

Fast response and easy interface usage also adapt flawlessly to all of your clinical needs, giving your assessments comfort and flexibility.











Fast and easy



21.5" LCD widescreen monitor



Zero-click automation



Tailored configuration



Large probe portfolio

- ✓ Battery
- ✓ Booting time less than 15 sec\*

# Large probe portfolio

Transducers are the core of Ultrasound technology. Integrating physics, electronics and geometrics in their design is **the greatest engineering challenge** of the Signal Processing Chain.

Transducers are the primary component of a Signal Processing Chain, the system that leads to the final diagnostic image. Although a great deal of time has been spent on the optimization of scan converters, post-processing algorithms, and sophisticated speckle-reduction technologies, **ultrasound transducers** remain a scanner's primary interface between patient and user.





The design, material, and manufacturing technology of transducers are the main determinants of an ultrasound system's image **quality**. Thanks to the innovation of gold standard ultrasound transducers, iQProbes offer state-of-the-art imaging.

- ✓ Active matrix composite material
- Single crystal
- ✓ Multiple adaptive layers
- ✓ Bi-con geometric lens
- appleprobe design
- Extensive use of applications with extended wideband convex, linear, phased array, volumetric, intraoperative and special transducer shapes.


\*from stand-by mode

## Clinical tools

### Zero-click automation


# Applications

#### ElaXto



ElaXto is a non-invasive method that supports the physician in assessing tissue elasticity. The differences in tissue responses are detected and visualized in real time.

#### Needle visibility



Enhanced and clear visualization of the needle during intervention procedures.

#### Stress echo



Complete Stress echo package with flexible and customizable protocols for imaging acquisition and review, also available with

#### XLight



Advanced algorithm to improve volumetric rendering

#### AutoNT



Automatic measurement of Nuchal Translucency (NT).




Automatic measurement of the Ejection Fraction (fully

#### QIMT



Automated real-time detection of Intima Media Thickness, including standard deviation and reliability index, based on RF signal analysis.

#### XStrain™



Global strain bullseye (17 segments) as a result of the 3 apical GLS outcomes.

## AutoEF



automated).

# Women's health



The convex and endocavity probes provide excellent image quality for women's health applications.

The 3D convex probe can also be used for standard examinations.

# Cardiovascular



MyLab™X6 is equipped with comprehensive cardiac and vascular configurations. It is a complete system for any cardiovascular ultrasound exam featuring customizable measurements and reporting.

# (1/1) General imaging



Esaote's new MyLab™X6 covers all clinical needs, from abdominal to endocrinological applications, to establish a diagnosis and provide the best possible therapy and follow-up.













Please visit us online for more information



Esaote S.p.A. – sole-shareholder company Via Enrico Melen 77, 16152 Genova, ITALY, Tel. +39 010 6547 1, Fax +39 010 6547 275, info@esaote.com